**Mesures**

***Measurements***

**Réf :  
282 009**

**Français – p 1**

**English – p 7**

**Version : 2104**

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***Power level indicator***

## 1 Généralités

### 1.1 But de l'appareil

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# 1 General

## 1.1 Purpose of the instrument

The power level indicator is used to make a qualitative evaluation of a power level reached.

It can thus be used to compare electrical powers available at the terminals of different alternators rotated by:

- a mechanical drive,
- falling water,
- an air jet.

It can thus:

- demonstrate the generation of electricity by an alternator,
- become familiar with parameters used in an «electricity power station», to optimise the available power.

## 1.2 Description

①	LEDs indicating a power level
②	Reset the memorized maximum power indication
③	Inclination prop to improve visibility
④	Potentiometer to obtain a power level indicator adjustment at full scale
⑤	Input voltage through Ø 4 mm safety bushings Maximum input voltages: 30 VAC and 15 VDC
⑥	Housing for 6LR61 9 V battery (provided)

## 1.3 Technical characteristics

### 1.3.1 Operating principle

The input voltage is compared with a reference voltage and displays the result by lighting up LEDs (from 1 to 11), proportionally to the ratio between the input voltage / reference voltage.

The reference voltage is corrected to zero in order to obtain a significant illumination.

The input signal, rectified by diodes and then amplified, is input into two special integrated circuits installed in cascade. This signal is amplified using a potentiometer (adjustable with a screwdriver) to obtain a full scale adjustment.

A variation of the input voltage causes LEDs to light up or go off according to a linear scale.

The last LED remains lit up to show the maximum level that was reached. A reset (RAZ) linearly reduces the maximum voltage reached to 0 V.

Therefore, the result is an analogue display device capable of displaying the variation of the input voltage.

The power level indicator is not a wattmeter. It is a simple instrument that qualitatively provides information about power concepts.

### 1.3.2 Usage conditions

The power level indicator is specifically used for experiments carried out with the associated equipment: reversible demonstration alternator (ref. 292 075), demonstration alternator model (ref. 292 077), water and air turbine kit (ref. 282 006), ecological lamp (ref. 332 016), hydroelectric power station (ref. 282 005), etc. It can also be used for other voltage sources less than 30 V AC and 15 V DC.

The advantage of the power level indicator is that it can make a comparison of **electrical powers available** at different alternator terminals or at alternators moved by different mechanical actions.

Thus, the comparison is only meaningful if the power level indicator is connected to the terminals of each alternator at no load (the current that passes through the power level indicator is of the order of a few micro-amperes). The electrical power available at the alternator terminals is thus proportional to the induced e.m.f.

In the case in which the power level indicator is used under load (for example incandescent lamp), the **electrical power actually output** to the circuit is displayed. This is less than the available power due to Joule effect losses in the circuit (load and resistance of the alternator winding).

#### In conclusion :

The power level indicator displays:

- at no load, the available power (potential).
- under load, a real power.

These two powers cannot be compared.

## 2 Usage

### 2.1 Full scale adjustment

The full scale adjustment allows to adapt the power level indicator to the tested equipment. It requires a small screwdriver.

Procedure to be followed :

- connect the equipment,
- realize a manipulation giving a maximal signal,
- when the signal displays on the LED, turn the potentiometer ④ using the screwdriver to adjust the full scale.

### 2.2 Example usage cases

#### 2.2.1 Reversible demonstration alternator (ref. 292 075)

#### 2.2.2 Demonstration alternator model (ref. 292 077)

#### 2.2.3 Air and water turbine kit (ref. 282 006)

## 2.2.4 Ecological lamp (ref. 332 016)

### Note:

Two experiments can be carried out:

- 1<sup>st</sup> experiment: with the incandescent lamp unscrewed from its base, rotate the crank handle at about 1 turn/second, providing a certain muscular effort (driving torque C1). The 9<sup>th</sup> LED lights up.
- 2<sup>nd</sup> experiment: with the incandescent lamp screwed into its base, turn the crank handle, once again at about 1 turn/second, providing a «different» muscular effort (driving torque C2). The 4<sup>th</sup> LED lights up

Obviously, the muscular effort supplied under load is greater than it was at no load:  $C2 > C1$ . And yet only four LEDs light up!

In the 1st experiment, the power level indicator provides information about the available power (potential), while in the second experiment it indicates the real supplied power.

The qualitative justification is that the muscular effort (mechanical energy) to be supplied actually depends on the electrical energy consumed. Thus if there is no lamp, there is no consumed electrical energy and therefore no mechanical energy to be supplied (except for the energy basically to overcome resistance and friction in the gear system).

This qualitative approach can be extended to include daily realities:

- Does a car use more fuel when its headlights are on?
- Can «electricity» be stored?

### 2.2.5 Hydroelectric power station (ref. 282 005)

### 2.2.6 Example use: Comparison of two powers

The output power levels from the reversible alternator (at no load) ref. 292 075 and the ecological lamp (at no load = incandescent lamp unscrewed), ref. 332 016 are evaluated, with the same sensitivity setting for the power level indicator.

It is observed that more LEDs light up with the ecological lamp than with the reversible alternator (no matter how fast it turns).

It is deduced that the ecological lamp alternator outputs a higher power than the reversible alternator.

#### Explanations:

The former generates a high current and it can light a LED and an E10 bulb up to 6 V / 3 W.

The latter generates sufficient current to light a LED but not enough to light up an E10 type bulb.

We could ask the question: why is the power greater in one case than in the other? We could then look for clues to the answer:

- coils,
- magnets,
- rotor rotation speed,
- etc.

### **3 After-Sales Service**

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries