**Méthode ESAO<sup>®</sup>**

***ESAO<sup>®</sup> method***

**Ref :  
453 040**

**Français – p 1**

**English – p 5**

**Version : 6006**

**Fermenteur  
*Fermenter***

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Example of the experimental device  
(probes not supplied)


## 1 Presentation

The fermenter is a cell culture chamber to study the metabolic pathways (respiration, fermentation, photosynthesis). It helps carry out long experiments.

It was specially designed to carry out Computer Assisted Experiments.

## 2 Use

- After calibration, place the probes in the holes using the pierced yellow corks,
- Plug the unused holes with the solid yellow corks,
- Wrap the top of the immersion heater in a strip of woollen cotton and insert it in the hole with the largest diameter. The role of this cotton seal is to prevent an increase in pressure in the fermentation chamber as well as to avoid its contamination from external pollution.

- Pass the aerator pipe through the centre hole of the lid and slide it to adjust the position of the diffuser.
- Position the centring disc at the bottom of the tank and place the tumbler in the centre

Place the medium in the glass chamber (yeast, chlorella, etc.); you require about 1.5 litres of the preparation.

- Install the 2 syringes on the valve and the pipe on the third channel of the valve. Pass the pipe through a passage of Ø 6 mm on the lid. One of the syringes may be used to inject a reagent and the other to take a sample of the preparation. Make sure to drain the pipe by flushing it several times.
- Close the chamber by putting on the lid attached with the PVC seal.

### 3 Composition


The fermenter consists of a chamber of about 2 litres with its “suitable lid”. It is delivered with:

- 4 solid corks for 12 mm holes
- 4 pierced corks for 12 mm holes
- 4 O-rings Ø 12 mm
- 2 O-rings Ø 6 mm
- 2 solid corks for 6 mm holes
- 2 pierced corks for 6 mm holes
- 1 cork to plug the opening for the immersion heater
- 1 immersion heater for the temperature control of up to 32 °C
- 1 150 l/hour aerator with diffuser
- 1 sampling device consisting of 2 syringes, a 3-way valve and a pipe
- 1 tumbler for stirring and a centring disc for the tumbler at the bottom of the tank
- 1 seal placed between the tank and the lid
- Some woollen cotton

### 4 Example of results

An example of results obtained during fermentation of yeast culture with the Sérénis software. Demonstration of the consumption of oxygen and the production of CO<sub>2</sub> and Ethanol.

## 5 Precautions for use

Precautions for use relating to the immersion heater: Read the instructions carefully. The immersion heater must be connected to the mains only when it is immersed up to a height of 15 cm. Operating it in air may result in permanent damage.

Do not wash in a dishwasher.

Do not sterilise it with spirit. Sterilise it with diluted Javel water.

## 6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
Rue Jacques Monod  
BP 1900  
27 019 EVREUX CEDEX FRANCE  
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

**Free service \* :**  
**+ 33 (0)2 32 29 40 50**

*\* Call cost not included*