ExAO - *Measure*

CAE - *Measure*

**Ref :
453 001**

Français – p 1

English – p 11

Version : 2104

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
1 Description


2 How it works


This probe works according to the principle of a Clark electrode. It measures oxygen amperometrically. The probe works **in a liquid or gaseous medium**. The electrode consists of a platinum cathode and a silver anode immersed in an electrolyte (50 g/L alkaline sodium phosphate [Na₃ PO₄] solution).

A potential difference of 800 mV is established between the cathode and the anode, and the oxygen present between the electrodes is reduced according to the following reactions:


- at the cathode:

- at the anode:

- the overall reaction:

The intensity of the resulting current is proportional to the concentration of oxygen in the electrolyte.

This current is amplified, and after calibration, the percentage or the mg/L of oxygen in the medium can be determined directly from it.

Under the proper conditions for use, the response of the probe is a linear function relative to the oxygen concentration in the medium. The intensity measured follows the law:

$$i = a \times [O_2] + b$$

i = intensity measured

a = calibration slope ("slope")

$[O_2]$ = oxygen concentration in the medium

b = ordinate at the origin ("zero").

3 Content

- Probe body with sensor head and protective tube
- Electrolyte Solution
- Solution "zero" for calibration
- A head replacement probe

4 Supplement

This probe can be used with all interfaces Jeulin. It requires the use of an adaptor or a sensor

- | | |
|-------------------------------------|----------------|
| - Oxygen sensor (Foxy) | Item no 482039 |
| - Oxygen sensor (Primo, Synéo, VTT) | Item no 472101 |
| - Oxygen adaptor (Visio) | Item no 482155 |

With the aid of the connector, connect the probe to the adapter and plug the adapter into the corresponding receptacle on the interface.

5 Commissioning

- Carefully remove the protective tube. To do this, do not turn anticlockwise. Rotate in the direction of clockwise or better pull gently.
- Unscrew the probe body
- Fill the probe electrolyte in three quarts
- Remove air bubbles included in the electrolyte

- Screw the probe head so that the electrolyte overflows. This ensure the absence of air
- Connect the probe to the sensor or adaptor
- Turn on the interface
- Wait 10 minutes (polarization time of the probe)
- Refer to manual of the sensor or adaptor used to calibrate the probe

6 During the experiment.

Note: the membrane of the probe head is very fragile, it should be handled carefully.

- Always remove the protective tube by turning it clockwise in order not to unscrew the probe head
- In the bioreactor, place the probe on one side and not at the center to avoid contact with the magnetic stirrer
- At the end of the experimentation, rinse the probe with distilled water to remove small particles that could clog the pores of the membrane

7 Storage

- **Between two experimentation:** It is not necessary to disassemble the equipment. Soak the probe as is in distilled water.
- **For a maximum period of one week:**
 - o Disconnect the probe adapter / sensor
 - o Place the probe in a beaker containing distilled water bottom
- **For longer than one week:**
 - o Empty the electrolyte to prevent oxidation of the probe or the evaporation of the electrolyte
 - o Rinse the electrodes and the probe head with distilled water
 - o Carefully wipe the electrodes
 - o Keep the probe body dry by protecting the electrodes using the protective tube
 - o Keep the probe head in a wet (plastic bag with a damp paper towel) to prevent drying of the membrane.
 - o The restart is done by following the procedure detailed in 5.
- **Once a year:**
 - o It is highly advisable to change the probe head once a year
 - o In case of black spots appear on the silver electrode, gently clean the oxidized parts with a soft cloth dampened with a solution of 0.1 M ammonia Rinse thoroughly with water distilled.

8 Malfunctions: causes and remedies

The measure remains perfectly constant:

- Check that the probe head is screwed: if the screw
- Verify that the membrane is not pierced: if there is a hole, quickly screwing the probe head a drop of solution. Electrolytic come out from that hole, then it is apparent on the membrane. If the membrane is pierced, replace the probe head.
- Check that the probe head contains many of the electrolyte solution: otherwise recover
- Verify the absence of bubbles in the probe head or in liquid medium in the probe head.
- Check that the electrode is not oxidized: this translates into a black deposit that can be cleaned with a dilute ammonia solution well.

If the measures appear attenuated or slow: in case of evaporation, the electrolyte solution is then supersaturated, crystals form and clog the pores of the membrane which results in malfunction: exchange surface between the middle and the reduced solution, resulting in slower exchange undervaluation variations. We must change the electrolyte solution

9 Technical specifications

Attention: Any substance diffusing through the membrane and having a reduction potential lower than 800 mV will react with the cathode and the anode; therefore, there may be interferences. Likewise, gases that react in oxygen (Cl_2 , SO_2 , H_2S) or those that affect the pH of the electrolyte (CO_2 , NH_3) will modify the values measured by the probe.

10 Spare parts

- Probe head..... Item no. 453 130
Factory-assembled, this head cannot be taken apart, thus assuring an optimum seal. The head fits directly on the body of the probe and is therefore easy to manipulate.
- 50 ml electrolyte bottle with dropper Item no. 106 084
- Zero solution
 - * ready to use, set of 2 x 25 ml ampoules..... Item no. 107 265
 - * pure anhydrous sodium sulfite Item no. 105 154
To be dissolved in water until saturation and then filtered

11 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50