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***Measurement and instrumentation***

**Ref :  
251 042**

**Français – p 1**

**English – p 3**

**Version : 8009**

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## 1 Technical specifications

- Power supply:
  - 1 6F22 type 9 V battery
  - or 1 12 V mains adaptor (not supplied)
- Ranges: -40 to +120 °C and 0 to 1000 °C
- Accuracy:
  - 2% of the value read +2 resolution units, excluding probe accuracy
- Resolution:
  - 0.1 °C on the -40 to +120 °C range
  - 1.0 °C on the 0 to 1000 °C range
- Probe:  
K type thermocouple (not supplied)
- Analog output:
  - -5 V to +5 V proportional to the temperature for the 2 ranges, available only when the device is powered by a 12 V mains adaptor. Accuracy: 1% of the displayed value
- Dimensions:
  - 100 x 100 x 40 mm

## 2 Heating time of the device

After powering on, the device requires about 7 minutes to attain its maximum accuracy. The "Cold" error is of at least 1 to 1.5 °C in relation to the actual value.

## 3 Replacing the battery

It is recommended that the battery be replaced when the "BAT" indicator flashes at the top left of the display. To open the battery housing, insert a small screwdriver in the notch at the edge of the flap and lever it up.

## 4 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**0825 563 563**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries