Microscopie

Microscopy

Réf :
571 210
571 211
571 212
571 261


Français – p 1

English – p 9

Version : 1111

Microscopes IRIO MS / ES / EPS / EPSA

Microscopes IRIO MS / ES / EPS / EPSA

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
0 825 563 563*

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Model IRIO EPS Reference 571 212

- A** Monocular head
- B** Prehension zone
- C** Stand
- D** Macrometric movement and stop
- E** Micrometric movement
- F** Indexing of the micrometric movement
- G** Compartment of the transformer unit
- H** Large field eyepiece x10 secured
- I** Turret objective holder
- J** Secured objectives
- K** Ergonomic stage
- L** Abbe's condenser numerical aperture 1.25
- M** Diaphragm and filter holder
- N** Cold integral lighting
- O** Adjustment of the lighting intensity
- P** Annotation zone

1 Introduction

IRIO: Another vision of microscopy

- Developed step by step with a group of teachers,
- Scalable: you can upgrade your equipment depending on your budget,
- Four patents: truly innovative technical solutions that are validated for the safety, durability and fast run through of your microscope.

Three versions that adapt to your environment and budget.

At any moment you can transform your IRIO-MS or IRIO-ES to IRIO-EPS.

2 Composition

I: included components, A: adjustable component

| COMPONENTS | | | Microscope | | |
|-------------------------|---------|--|------------|---------|-----------------------|
| TYPE | REF | DESCRIPTIF | IRIO MS | IRIO ES | IRIO EPS IRIO EPSA |
| Stand | | | | | |
| Stand IRIO | - | - Turret 4 objectives - fast and fine focussing | I | I | I |
| Mirror | - | | I | - | - |
| Integral lighting | 574 241 | | A | I | I |
| Ergonomic stage | 574 242 | | A | A | I |
| Observation head | | | | | |
| Monocular | - | - Inclined at 45° and pivoting through 360° | I | I | I |
| Eyepieces | | | | | |
| Large field x10 secured | - | - WF x10 /18 | I | I | I |
| Large field x15 secured | 574 244 | - WF x15 / | A | A | A |
| Huyghens x5 secured | 574 243 | - H x5 / | A | A | A |
| Objectives | | | | | |
| Achrom. x4 secured | - | - x4 O.N. 0,10 | I | I | I |
| Achrom. x10 secured | - | - x10 O.N. 0 ,25 | I | I | I |
| Achrom. x40 secured | - | - x40 O.N. 0,65 retractable | I | I | I |
| Achrom. x60 secured | 574 246 | - x60 O.N. 0,85 retractable | A | A | A |
| Achrom. x100 secured | 574 247 | - x100 O.N. 1,25 retractable, immersion | A | A | A |
| Complements | | | | | |
| Adjustment key | 574 245 | Assembly/disassembly/stop adjustment | I | I | I |
| Integral lighting | 574 241 | Replacement of mirror | A | I | I |
| Ergonomic stage | 574 242 | Replacement of standard stage | A | A | I |

3 Activation

1. Remove the microscope from the packing.
2. The microscope that you receive is pre-assembled and ready to use. The eyepiece H and the objectives J are in position and tightened (secured).
3. For models with integral lighting (Ref. 571 211 and 571 212) connect the cord with the 2P + T to a 220 V main plug.
4. Keep the key, dust cover and polarisation kit that are delivered with your microscope.

4 Utilisation and adjustments

4.1 Adjustment key

This tool, which is provided with each IRIO microscope, allows the user to carry out the following actions:

- Disassembly/locking of the eyepiece **H** or objectives **J**.
- Opening of the macrometric movement **D** and adjustment of the stop.
- Disassembly of the mirror unit or lighting **N**.

In case of loss, this key can be ordered under the reference 574 241.

4.2 Stand

- Ergonomic stand **C**: stable pedestal with smooth and comfortable shape, easy to store.
- Excellent prehension by means of a robust beam **B** that is used as a handle.

- Annotation zone **P**: the identification of each microscope is clearly visible.

4.3 Observation head

- Monocular head **A** inclined at 45° and rotary over 360° suitable for observation by two users.
- Eyepiece x10 secured large field **H**: benefits from fall-prevention, dustproof and patented antitheft devices. The disassembly remains possible with the provided key.

The eyepiece holder is standard and can hold eyepieces Ø 23 mm (DIN norm). This eyepiece is compatible with the cameras and contains a compartment to place the analyser (provided) in polarised light.

Optional: eyepieces x5 Huyghens (Ref. 574 243) and x15 large field (Ref. 574 244).


Change the eyepiece:

Place the key in the flat plate zones.

Unscrew the eyepiece while firmly maintaining the key with your other hand. You can use a latex glove or a cut piece of sheet rubber.

4.4 Objectives

Can only be disassembled with the provided key. These objectives will be protected from deterioration, dust and theft. Each microscope is delivered with achromatic objectives x4, x10 and x40 (retractable) on turret I 4 places. Magnifying power range x40 to x400 with the provided eyepiece x10.

Optional: objectives x60 (retractable) Ref. 574 246 and x100 (retractable, immersion) Ref. 574 247.

Change the objective:

Place the key on the ring openings opposite the flat plate zones on the objective. Unscrew the objective.

Flat plate zone to place the key

4.5 Focussing - Stop


The helical type macrometric movement **D** replaces the toothed bar and guarantees your comfort: no sustenance and wear. **WITHOUT A TOOTHED BAR** you are now freed from the sustenance and maintenance of the macrometric movement.

Recommended position of use:

Adjustment of the stop:

For the protection of the preparations: a stop **D** can be obtained from the teacher and can be adjusted with the provided key.

Introduce the key on the edge of the macrometric movement **(1)** and open **(2)**. Place the stage to the maximum desired height (objective of greatest magnifying power). Place the stop screw in contact with the axis of the spiral and tighten with the key **(3)**.

Important: for this action, position the micrometric movement as shown in the picture below.

Important: Never change the adjustment of central screw: it is the adjustment of movements' hardness which is factory-adjusted, you may damage your material. For any adjustment of movement's hardness, contact the JEULIN technical support.

The micrometric movement **E** is manipulated in a typical manner. An indexing **F** of the min/max excursion provides additional safety: when the left mark is positioned in the middle of the two right marks the micrometric movement is at half-lift.

4.6 Stage

| | Microscope MS Ref. 571 210 | Microscope ES Ref. 571 211 | Microscope EPS Ref. 571 212 |
|--------------|-------------------------------|-------------------------------|--------------------------------|
| Stage | 120 x 140 mm with bases | 120 x 140 mm with bases | 120 x 140 mm ergonomic |

Depending on the IRIO model of the microscope, you are provided with either a rectangular stage equipped with bases (microscopes Ref. 571 210 and 571 211), or an ergonomic stage (microscope Ref. 571 212).

Intuitively practical and precise, the ergonomic stage incorporates a plate holder for histology plates and rock thin sections and commands with the buttons of orthogonal movement. It is designed so that the pupils can quickly run through it.

Position of use:

Push the plate holder towards the left. Place the plate in the plate holder and release gently.

Move the plate, either roughly with the stage, or finely by means of the two buttons.

Ergonomic stage K (Reference 574 242)
 The use of the bases is difficult when the magnifying power is high, therefore you can upgrade your microscope by incorporating the ergonomic stage. Practical and precise, it guarantees easy use by pupils. The stage is delivered with the necessary fixtures and a detailed manual.

Unscrew the four retaining screws of the standard stage, remove the wires and place the ergonomic stage without spacers (see manual delivered with ref. 574 242).

4.7 Condenser, diaphragm, filter holder and polarisation

Depending on the IRIO model of the microscope, you are either provided with a rectangular stage equipped with an iris diaphragm and a filter holder **M** (Ref. 571 210 and 571 211), or an iris diaphragm, a filter holder and an Abbe's condenser **L** (Ref. 571 212).

- Iris diaphragm continuously adjustable.
- Filter holder: it allows for the positioning of the polarizer (provided) when thin sections are studied in polarised light. A blue filter can also be placed in the filter holder but the integral lighting **O** of the IRIO microscopes does not require this.
- Abbe's condenser, numerical aperture 1,25 (only for the reference 571 212), incorporated and adjustable, height can be set by means of a handle (helical movement).

Polarisation kit: consists of a polarizer that must be placed in the filter holder and an analyser that must be clipped on top of the eyepiece (push on the small red pearl).

Polarizer

Analyser

Analyser in place

4.8 Lighting

| | Microscope MS Ref. 571 210 | Microscope ES Ref. 571 211 | Microscope EPS Ref. 571 212 |
|----------|-------------------------------|-------------------------------|--------------------------------|
| Lighting | With mirror | With LEDs | With LEDs |

Depending on the IRIO model of the microscope, you are either provided with lighting by a mirror (microscope Ref. 571 210), or with lighting with incorporated leds (microscopes Refs. 571 211 and 571 212).

To use the mirror, the ideal position of the add-on lamp is sideways, therefore do not place your hand between the lamp and the mirror (see picture opposite).

- Integral lighting **N**: you benefit from good quality, cold, durable lighting and its intensity is adjustable with the potentiometer **O**. The emission spectrum of the LEDs (LED indicators) does not require the addition a blue filter.

It is recommended to begin the observation always with the luminous intensity set to its minimum (button stopped towards the left).

Lighting block (Reference 574 241)

This low voltage built-in lamp allows for the transformation of your mirror microscope IRIO-MS to a lighting model without any technical assistance. This easy to assemble kit contains:

- A lighting unit with the LED bulbs, the adjustment potentiometer, and an outlet,
- A power supply with jack connection to the lighting unit,
- The fixture screws and a protection mask,
- A detailed manual with pictures.

Tools required: the IRIO adjustment key and two screwdrivers.

Remove the mirror unit and clip the lighting unit at the same location. Turn the microscope down and install the power supply unit (see the manual provided with the ref. 574 241).

- Power supply: built-in transformer unit **G** 12 V / 160 mA / 2 VA located in the stand **C**: no encumbrance of the low voltage wall socket.


5 How to settle (adjust) the hardness of the movement for microscopes IRIO?

- Introduce the key on the slice (edge) of the macrometric movement and open it.
- Fall saw her(it) central **(1)** clockwise.

Attention: this operation is only made if saw her(it) central is steel (silvery) and not in brass (gilt). If saw her(it) is in brass, it is necessary to turn(return) the microscope with our SAV (AFTER-SALES SERVICE).

- Around this lived central, there is a copper-colored slice **(2)** provided with 2 small openings. Arm yourselves with a small screwdriver and with a hammer and tighten slowly this slice clockwise.

- It is then necessary to re-screw saw her(it) central counterclockwise.

6 Maintenance and precautions of use

- The microscopes must be placed far from water projections and preferably in a room with atmosphere moisture content lower than 70 %.
- Do not store the microscopes in a cupboard where reagents that can emanate aggressive vapours are stored.
- To avoid dust to deposit, it is recommended to leave the eyepieces in their tubes and to cover the device with its plastic dust cover after each use.
- Dust that may be left on the optical parts must be blown or removed by means of a dust collector (Ref. 574 198).
- Avoid touching the surface of the lenses with your fingers.
- In order to clean finger marks, use a soft lint-free cloth or a cleaning paper (Ref. 803 092).
- After observation with immersion objectives, the oil that is stuck to the front lens of the objective should be gently wiped off with a soft cloth. Never use alcohol, but use a specific cleanser for optics (Ref. 102 146).

7 Technical characteristics

For the models with lighting (Ref. 571 211 and 571 212).

- Supply voltage: 230 V - 50 Hz
- Power: 5 VA
- Class 2 device
- Environmental conditions:
 - Use exclusively indoors
 - Ambient temperature between 10 °C and 40 °C
 - Relative moisture content lower than 70 %
 - Fluctuation of the power supply network not being over ± 10 %
 - Installation category II
 - Degree 2 of pollution
 - Altitude of use < 2000 m
- The power supply plug is used as a disconnecting switch to cut the energy supply of the device.

Comment: if the device is used differently than specified below (section « utilisation and adjustments»), the protection provided by the device may be compromised.

8 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques-Monod
BP 1900
27019 EVREUX CEDEX France
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries