Didactique

Kit de chimie

Didactic

Chemistry kit

Ref :
105 009
106 044

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1 Description

1.1 General

Constituted of a set of pre-dosed solutions, the “ion characterization” test kit allows substantially reducing the time of preparing the TP for ion characterizing tests by streamlining the equipment used to make the dilutions and the glassware present in the laboratory. The solutions supplied are designed for correct visualising of a precipitate in a solution. During the experiments students observe that they are working in two distinct media.

The quantities of ionic compounds present in the dropper bottles are reduced to a minimum to diminish risks linked to these products in the laboratory. By preferring plastic packaging, the risk of breakage, cuts and the cost of recycling the packaging for chemical products are reduced to a minimum.

1.2 Description and composition

1.2.1 Secondary school ion characterization test kit (ref. 106 044)

The ion characterization test comprises:

5 plastic 30 mL dropper bottles containing:

- A solution of copper sulphate
- A solution of iron II sulphate
- A solution of iron III sulphate
- A solution of zinc nitrate
- A solution of sodium chloride

1 yellow glass 30 mL dropper bottle containing:

- A solution of silver nitrate to identify the anion.

1 plastic 60 mL dropper bottle containing:

- A sodium hydroxide solution for identifying the cations.

This allows performing 6 tests ten times for the following ions: Ag^+ , Cu^{2+} , Fe^{2+} , Fe^{3+} , Zn^{2+} and Cl^- .

1.2.2 Grammar school ion characterization test kit (ref. 105 009)

The ion characterizing test kit comprises:

9 plastic 30 mL dropper bottles containing:

- A solution of copper sulphate
- A solution of iron II sulphate
- A solution of iron III sulphate
- A solution of zinc nitrate
- A solution of sodium chloride

- A solution of calcium chloride
- A solution of ammonium oxalate to identify the calcium ions.
- A solution of barium chloride to identify the sulphate ions.
-

1 yellow glass 30 mL dropper bottle containing:

- A solution of silver nitrate to identify the anion.

1 plastic 60mL dropper bottle containing:

- A solution of sodium hydroxide to identify the cations.

This allows testing ten times for the following ions: Ag^+ , Cu^{2+} , Fe^{2+} , Fe^{3+} , Zn^{2+} , Cl^- , SO_4^{2-} and Ca^{2+} .

The reagents used to create the reactions are sodium hydroxide (Na^+ , HO^-)aq, silver nitrate (Ag^+ , NO_3^-)aq., barium chloride (Ba , Cl^-) and ammonium oxalate monohydrate ($(\text{NH}_4^+)_2\text{C}_2\text{O}_4^{2-}$)aq.

2 Additional equipment

To conduct all the experiments the following equipment is required:

Designation	Secondary school	Grammar school
Basic rack for 6 tubes	1	2
Test tube	6	8
Pair of protective gloves in vinyl or latex (1)	1	1
Pair of safety glasses (1)	1	1

Table 1

(1) Individual protection equipment recommended although not mandatory

3 Usage and safety precautions

3.1 Warning

The Jeulin Company cannot be held responsible for an accident that occurs after the iron protection and corrosion kit has been used under conditions other than those described in these instructions.

Similarly, the Jeulin Company cannot be held responsible for an accident that occurs because the safety instructions given in these instructions are not respected.

3.2 Labelling and safety datasheet

Flasks containing the reagent are labelled individually in accordance with the regulations in force.

Before performing any manipulation, get pupils to read the label so that they can be familiar with the risks and protection measures to be applied.

If any additional information is required, the safety datasheet for the solution can be downloaded from the www.jeulin.fr site

Product	FDS
Ions specific test kit – silver nitrate solution	106 044-A
Ions specific test kit – copper sulphate solution	106 044-B
Ions specific test kit – Iron II sulfate solution	106 044-C
Ions specific test kit – Iron III sulfate solution	106 044-D
Ions specific test kit - zinc nitrate	106 044 E
Ions specific test kit – ammonium oxalate solution	106 044-F
Ions specific test kit – barium chloride solution	106 044-G
Ions specific test kit – calcium chloride solution	106 044-H
Ions specific test kit – sodium hydroxide solution	106 044-I

FDS = Safety Datasheet

3.3 Individual protection

The iron protection and corrosion kit was formulated to minimise risks related to use of the product. Although not compulsory, we recommend that safety goggles shall be worn in addition to the lab coat during handling. It is also recommended that gloves should be worn.

3.4 Collective protection

No collective protection equipment is required during storage (safety cabinet, etc.) or during use (hood, etc.), because reagents do not contain any flammable volatile or toxic solvents.


3.5 Waste

The iron protection and corrosion kit does not contain any components that create any particular danger for the environment. However, these products must be considered as being chemical laboratory waste and treated as such, in the context of a citizen environmental protection policy. Reagent quantities used in this kit are kept low in order to minimise waste volumes and disposal costs.


4 Procedure

4.1 Principle of ionic precipitation

The ion characterization test kit allows revealing an anion or a cation depending on the reagent used. During the reaction an anion of a chosen solution will react with a cation in the reagent and vice versa. The anion and the cation which do not react are called spectators. These spectator ions constitute the liquid medium which is insoluble vis-à-vis the precipitate formed. The principle of precipitation is based on mixing two saline solutions which are miscible. Contact between the positive and negative charges contained in the mixture of saline solutions gives rise to a solid non-miscible precipitate. The reaction equation follows the model.

In practice the experiment can be defined as follows:

A saline solution with the ion to be identified and a reagent is poured into the first test tube on the left. The reagent contains the ion with the reverse charge to the ion to be identified. During mixing of the two solutions, an insoluble precipitate is formed in the aqueous medium.

4.2 Preparation of the operation

Assemble at a single work station the recommended accessories (paragraph 2 and photo 1) and a bottle of each ionic solution.

Photo 1

4.3 Preparation of the characterization test

For the first use, cut the safety cover on the bottle (photo 2), after opening, pour into the test tube a volume which does not exceed 3 cm of the total capacity of the tube. To measure a constant volume for each experiment, a small graduated flask can be used.

Open in the same way if necessary, the bottle containing the reagent. Pour the equivalent of 1cm of the reagent into the tube (photo 3).

Photo 2

Photo 3

4.4 Observation of the precipitate

During mixing, the precipitate corresponding to the test is not yet formed. Mix the reaction medium and observe the precipitate formed (photo 4). The precipitate does not form immediately, to allow the student to observe the two media:

1st medium: A liquid medium incorporating the spectator ions.

2nd medium: A solid medium: the precipitate.


Photo 4

5 Cleaning of equipment

Immediately clean the test tubes after the operations to avoid any staining.

6 After-Sales Service

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries