Optique

Optics

**Réf :
201 016**

Français – p 1

English – p 4

Version : 8110

**Laser 2 mW
*2 mW Laser***

1. Description - Principe

1.1 But de l'appareil

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1. Description - Principle

1.1 Purpose of the instrument

A laser (abbreviation for Light Amplified Stimulated Emission of Radiation) is the ideal source for easily and immediately obtaining:

- A narrow and very powerful beam for geometric optical studies
- A beam with very high spatial and temporal coherence for enlightening experiments in physical optics using simple accessories.

The 2 mW laser is powerful and facilitates visibility of the experiments even further. Just connect to 220 V mains and you have the laser beam.

It is compliant with the safety rules defined by the circular of the Ministry of National Education dated 12/10/99.

WARNING:

Never look directly at a laser beam, or those reflected or refracted, or not sufficiently scattered.


The beam presents no risks for other parts of the body.

1.2 Tube technology

Coaxial design was chosen because of the very high stability of both power and direction that it gives the laser. In addition, the output mirror has a corrective lens on its front surface that is adapted to the geometry of the capillary and thus causing only minimal divergence. The cathode of the tubes is oxidized by a special process that increases their lifetime and enables problem-free storage. All the tubes are cured at 350°C to eliminate all residual gas impurities.

The ratio of the He-Ne mixture and total pressure have been optimized as a function of capillary diameter in order to obtain the highest laser power possible with the lowest noise.

- Size: 68 x 90 x 370 mm
- Weight: 1.4 kg

1.3 Technical characteristics

Power supply	220 V \pm 10%
Wavelength	632.8 nm
Minimum output power	2 mW
Maximum output power	5 mW
Beam diameter	Φ 0.75 mm
Beam divergence	\leq 1.2 mrad (10 mm at 10 m)
Polarization	Not defined
Priming voltage	\leq 7 kV
Operating voltage range	1200 to 1400 V
Operating current	5 mA
Break-induced current	\leq 3.5 mA
Minimum required ballast resistance	60 k Ω
Longitudinal modal space c/2L	614 MHz
Power stability (after heating for 20 min)	$\approx \pm$ 2.5%
Power when turned on	\geq 80%
Operating temperature	-25 to +80°C

1.4 Maintenance

This laser requires no special maintenance and can remain for several months without being used with no loss in operating characteristics.

Comment: in spite of the three anti-reflection layers on the output mirror, there is often a second main beam. This beam is less concentrated and much dimmer than the main beam (more than 10,000 times less). It thus does not interfere with optical experiments.

2. Experiments

2.1 Geometric optics

The qualities of a very fine and powerful 2 mW laser beam are used.

A light beam is obtained with a diffraction grating.

A broad beam can be obtained with the eyepiece of a microscope (a very low magnification objective or a lens with a 5 or 10 cm focal length).

Beams are visualized in air with smoke or in water by adding several drops of milk.

2.2 Physical optics

The accessories for laser optics are used and the experiments involve diffractions and interferences. Refer to the instructions accompanying these accessories.

3. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux