Mesure

Measurement

**Ref :
211 026**

Français – p 1

English – p 3

Version : 8009

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Large display lightmeter

1 Généralités

1.1 But de l'appareil

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1 Generality

1.1 Purpose of the instrument

The large display digital lightmeter is useful for making simple and fast measurements of real illumination.

1.2 Description

The lightmeter is composed of three parts :

1. A sensor formed of a photodiode and provided with a shutter.
2. A box containing the controls (On/Off switch, measurement range selector, memorise the current measurement) and display on an LCD screen.
3. A connecting cable.

1.3 Electrical specifications

- Three measurement ranges: from 1 to 1 999 lux; from 2 000 to 19 990 lux and from 20 000 to 50 000 lux.
- Memorise data on the display.
- Read on a 50 x 28 mm LCD screen
- Battery condition indicator.
- Display if the measurement range is exceeded.
- Photosensitive sensor composed of a silicon photodiode. The sensor provided with a removable shutter is connected to the control and read box through a cable that can be extended up to 1 m for easy use.
- Cosine correction (indicates precise values regardless of the angle of incidence of light on the cell).
- Switch: - « hold » function to display the most recently measured value
- manual range selector.

- Dimensions: - box: 118 x 70 x 29 mm
- sensor: 100 x 58 x 25 mm.
- Total mass: 180 g.
- Battery: 6LR61 9V (delivered).
- Sampling frequency: 2 measurements per second.
- Response time: 1 s.
- Precision: $\pm 4\%$ on the 10 – 20 000 lux range.
- Deviation due to temperature: 0.1% / °C.
- Fidelity: $\pm 2\%$.
- Set delivered in a rigid plastic case.
- Standard ISO 9001:2000.

Measurement range (Lux)	Resolution	Uncertainty
0 – 1999	1 Lux	$\pm (5\% + 2d)$
2000 – 19990	10 Lux	$\pm (5\% + 2d)$
20000 – 50000	100 Lux	$\pm (5\% + 2d)$

2 Use

2.1 Reminder about lux

The lux (symbol lx) is the unit of light illumination. The lux corresponds to a light flux of 1 lumen (symbol lm) per square meter: 1 lux = 1 lm.m².

2.2 Example use in chemistry

The lightmeter can be used for studying chemical dynamics when a product appears or a coloured reagent disappears. According to the Beer-Lambert law, light intensity transmitted by the reaction medium depends on the concentration of the coloured species. Thus, the influence of the concentration of a reagent, temperature or a catalyst can be demonstrated.

2.3 Measurement examples

Place	Applications	Measurements (Lux)
School	Laboratory, CDI	1500 – 750
	Classroom	750 – 300
	Inside a gymnasium	300 – 75

3 Replacing the battery

The battery must be replaced when the  symbol appears on the LCD screen.

- Remove the cover from the back of the instrument.
- Replace the battery. (6LR61 9V)

4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries