**Méthode ESAO<sup>®</sup>**

***ESAO<sup>®</sup> method***

**Ref :  
251 016**

**Français – p 1**

**English – p 4**

**Version : 7006**

**Sismomètre Initio SA**  
**Initio SA seismometer**

## 1 Introduction

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**NOTES**


## 1 Introduction

The seismometer allows displaying and recording a wave generated by a mechanical shock.

- You can **display the signals in real time** in the form of an LED which lights up increasingly brightly, depending on the intensity of the signal.
- You can **visualise and save the signal** in the form of an LED which lights up more brightly, depending on intensity of the signal.
- You can **transfer the raw signal of the wave to a computer assisted experiment interface** via an analogue  $\pm 5$  V output.

## 2 Description

The Initio SA seismometer is composed of a piezometric sensor linked to an Initio unit powered by a 9 V battery (not supplied).

- (2) The seismometer provides 2 modes:
- Upper position: real time display of signal.
  - Lower position: save signal.

## 3 Additional equipment

For a simple demonstration of the impact of distance on a wave, it is possible to use a seismic scale model (Ref. 507 036).

## 4 Start-up and first use

- Insert a battery in the housing of the casing located in the upper section of the seismometer,
- Start by using the on/off button located under the seismometer casing.
- If you want to use 2 seismometers for a quantitative approach to the phenomenon:
  - place the 2 seismometers side by side,
  - cause a shock which is equidistant from the 2 sensors,
  - adjust the response using the adjustment button on the casing.
- Select the use mode you want (recording of the wave signal or real time display of the wave signal).
- Place the red sensor flat on the surface where the shock is to be made.
- Cause a shock.

Since the sensor is extremely sensitive, a small shock (pencil falling from a height of 10cm or breaking a piece of plastic on the seismic scale model for example) is sufficient to obtain correct display of the signal.

### **Example of experimental set-up with the seismic scale model:**

## 5 Precautions for use

- Remember to switch off the seismometer after use checking the ON diode has gone out.
- Do not tap too hard, this could cause rebounding of the sensor.
- If possible, attach the sensor to the support which must itself be stable to avoid rebounds.

## 6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

**Free service \* :**  
**+ 33 (0)2 32 29 40 50**

*\* Call cost not included*