Acoustique

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Acoustics

Speed of sound

**Ref :
453 127**

Clap sonore Synchro

Synchronised clap sound

Français – p 1

English – p 3

Version : 7005

1 Description

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
1 Description

The apparatus comprises 2 metal rods of approximately 25 cm, with at the end, a female banana jack to connect a single or safety cable and a sleeve for insulation and grip.

The impact of the 2 rods striking one against the other generates a sound, and closes the circuit resulting in the triggering of data acquisition by your Ex.A.O system.

2 Composition

- 2 rods

3 Operation

Link the 2 rods to the synchronisation device on your measuring instrument or your interface using the 2 banana - banana cables.

- For an ESAO 4 interface, use the 2 red and black sockets located on the top of the interface,
- For a PRIMO or VTT console, use the 2 red and black sockets located on the side of the interface.

4 Examples of experiments: calculating the speed of sound in Ex.A.O.

The purpose is to measure the speed of the sound caused by striking the 2 rods together. This speed is 340 m/s at 15°C.

The impact of the temperature is not negligible; the speed varies according to the formula:

$$c_{\text{air}} = (331.5 + 0.6 \cdot \theta) \text{ m/s where } \theta \text{ is the temperature in } ^\circ\text{C}.$$

4.1 Equipment

- an interface with an external synchronise input
- generalist software
- a sound metre sensor
- 2 banana – banana cables of adequate length (minimum 1 metre).

4.2 Assembly

- Link the interface to the computer.
- Start your software.
- Connect your sound sensor to your console.
- Slide the icon of the sound sensor to the ordinate (acquisition can be in Decibels or signals).
- Set the sensor to zero.
- Slide the alarm clock icon corresponding to the time, to the abscissa.
- Set the time 10 ms and select external synchronisation.
- List the synchronised clap sound device to the synchronise input using a cable of adequate length (at least 1 metre).
- Start acquisition by clicking on the green light.
- Strike the two rods one sharply against each other (this triggers acquisition).

4.3 Results: acquisition in dB

Example with a distance of 1 metre at 25 °C
The propagation time is measured as 2.837 ms.
The distance between the sensor and the shock emission zone is 1 metre,
hence speed is 352.5 m/s.
The theoretical value is 346.5 m/s, that is a deviation of 1.7%.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*