# **Observation**

Microscopes

# ***Observation***

*Microscopes*

Réf :  
571 227  
571 271


Français – p 1

English – p 4

Version : 1105

## **Microscopes à LED trinoculaire**

## ***LED trinocular microscopes***

### Légende :

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
**Caption :**

- |                                    |                                         |
|------------------------------------|-----------------------------------------|
| <b>A</b> Wide field eyepieces      | <b>H</b> Inversed quadruple nosepiece   |
| <b>B</b> Achromatic objectives     | <b>I</b> Double layers mechanical stage |
| <b>C</b> Square stage              | <b>J</b> Coarse focus                   |
| <b>D</b> Adjustable Abbe condenser | <b>K</b> Stage adjustment               |
| <b>E</b> Illumination              | <b>L</b> Fine focus                     |
| <b>F</b> Trinocular head           | <b>M</b> Dimmer                         |
| <b>G</b> Video head                |                                         |

## 1. Technical specifications

<b>Head</b>	Trinocular, 30° inclined, 360° rotatable, non detachable	
<b>Eyepieces</b>	Wide field x10, fixed	
<b>Interpupillar distance</b>	55 – 75 mm	
<b>Dioptric adjustment</b>	1 tube	
<b>Achromatic objectives</b>	Inversed quadruple nosepiece	
	Réf 571 227	Réf 571 271
	x4 O.N. 0,1 ; x10 O.N. 0,25 ; x40 O.N. 0,65 retractable	x4 O.N. 0,1 ; x10 O.N. 0,25 ; x40 O.N. 0,65 retractable ; x60 O.N. 0,85 retractable
<b>Magnification</b>	x40 - x400	x40 - x600
<b>Stage</b>	Double layers mechanical stage 125x125mm	
<b>Focusing</b>	Separate coarse and fine focusing adjustment	
<b>Blocking lever of the stage</b>	Yes	
<b>Condenser</b>	NA 1,25 adjustable Abbe condenser with iris diaphragm and filter holder	
<b>Accessories</b>	Blue filter and dust cover	
<b>Illumination</b>	LED with separate adjustable brightness	
<b>Power supply</b>	Rechargeable batteries and separate 5,5V-300mA transformer	
<b>Weight</b>	Around 5 kgs	

## 2. Set-up

1. Remove the microscope from the packaging, and the various protections items. Save all packaging for long-term storage or return shipment.
2. Screw the video tube.
3. Screw the objectives 4x, 10x, 40x and/or 60x onto the nosepiece.
4. Use and test the mechanics of your microscope to be confident with it.
5. Plug in your microscope with the transformer.

**Note :** 1) Your microscope had to be plug in an earth connection.  
 2) Make sure that the supply voltage is adequate for the nominal voltage of the microscope.

The batteries have an initial charge time of 5 to 6 hours. The standard charge time is of 4 hours for an autonomy of 4 hours.

## 3. Use and adjustment

- When moving the microscope, pick it up only with the inclined arm of the body and place your other hand under the base.
- Use the switch behind the base of the microscope to turn it on, turn the dimmer located on the base to obtain the desired light intensity.
- Place the prepared slide on the stage, cover slip towards the objective lens.
- The 360 ° rotatable head and the 30° inclined tube are ergonomically designed for prolonged use.
- The microscope is focused vertically by displacing the stage with a coarse adjustment knob for rapid movements and a separate micrometric screw for fine adjustments.

- To prevent the objective touch the prepared slide when focusing, turn the coarse focus knob until the prepared slide should be closed to 3 mm from the objective. Turn slowly the coarse focus adjustment knob until a clear image is obtained; then use the fine focus knob to make the prepared slide image as clear as possible.
  - A light condenser is found under the stage. The brightness may be modified when the diaphragm aperture is adjusted.
- Note : The diaphragm aperture is not designed to adjust the brightness ; the brightness is adjusted using the brightness adjustment knob.**
- Turn the nosepiece to change the magnification, and make sure (by hearing a "click") the objective is in the light way.

## 4. Maintenance

- Protection from dust : After use, place the cover over the apparatus. If there are dust deposits on the lenses, blow them away with a syringe or dust remover (item no. 574 198) or use a martin hair brush ; the lenses are fragile and can be scratched by dust.
- Fingerprints : Do not touch the lenses. If, however, they do have fingerprints on them, remove the fingerprints with a soft cloth or lens paper (item no. 803 092).
  - Humidity : Humidity, especially in combination with temperature, is particularly harmful for assembled optics. It is suitable to store the microscopes in a dry place, and under especially exposed climate conditions. The microscopes should be stored in their polystyrene box with a desiccation bag.
  - Disassembly : A microscope is a finely-tooled precision instrument. If it is necessary to take it apart (for maintenance), the instrument should be returned to JEULIN or a qualified repair shop. The base plate that seals the microscope may only be opened by an electrician because the body must remain protected according to electrical safety standards. Do not take the microscope apart unless it is unplugged.

The shelf life of the LED is of 10,000 hours. In case of non-function of the lighting and after checking that the batteries are charged, send the microscope back to Jeulin.

## 5. After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries