Didactique

Kit de chimie

Didactic

Chemistry kit

Réf :
105 008

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This kit has been designed to be used by students.
Manipulation simple and fast, are achievable by all participants.
Without risk, there is no need to use a hood, only the usual protective equipment are recommended.

1 Description

1.1 General

Thanks to its innovative and original packaging, the Nylon 6-6 synthesis kit in break-open phial means the teacher can substantially reduce preparation time, by minimizing the quantity of glassware required. In addition, the break-open phial are pre-dosed and ready for use. So the reagents can be distributed individually to the students without any need to measure their volume.


The quantities of reagents in the break-open phial have been reduced to a minimum to limit risks linked to handling of the products. The formulation of the nylon 6-6 preparation kit in break-open phial is organic solvent free and does not release volatile organic compounds (VOC) so that it may be used in a classroom without using a hood.

Made of plastic, the break-open phial can be opened by students without any danger or risk of cuts.

1.2 Composition and description

The preparation of nylon 6-6 kit in break-open phial is composed of a set of 10 phials of organic phase containing 10 mL of adipoyl dichloride solution and 10 phials containing 10 mL of 1,6-diaminohexane aqueous solution.

This kit allows performing 10 syntheses.

Picture 1

- | | |
|--|--|
| 1 – Scored head of ampoules | 3 – Phial containing the organic phase |
| 2 – Label in compliance with regulations | 4 – Phial containing the aqueous phase |

Table 1

2 Additional material

The synthesis of nylon with the Nylon 6-6 preparation kit in break-open phial requires the following accessories:

Name	Quantity
Small 50 or 100 mL beaker or glass pot	3
Glass agitator or rod	1
Spatula	1
Sheet of absorbent paper	1
Piece of iron wire or paper clip	1
Pair of protective vinyl or latex gloves	1
Pair of safety goggles	1
Dishwashing detergent	1

Table 2

Individual protection equipment recommended even though not compulsory

3 Precautions for use and safety

3.1 Preservation

Caution: Products sensitive to light, keep refrigerated and use within 12 months after receipt.

3.2 Warning

JEULIN cannot be held liable in the case of accidents occurring during use of the Nylon 6-6 synthesis kit under conditions other than those set out in this notice.

Also, JEULIN cannot be held liable in the case of accidents which occur on account of non-compliance with the safety instructions described herein.

3.3 Labeling and Safety Data Sheet

The break-open phial containing the reagents are individually labeled in compliance with the regulations in force.

Before any operation, have the students read the label so they become aware of the risks and the protective measures to be adopted.

For additional information, the material Safety Data Sheets for the two reagents can be downloaded from the site www.jeulin.fr.

Product	FDS Reference
Nylon 6-6 preparation kit in phial – organic phase	105008-A
Nylon 6-6 preparation kit in phial – aqueous phase	105008-B

Table 3

3.4 Individual protection

The nylon 6-6 preparation kit in break-open phial has been formulated to minimize the risks associated with the use of the product. However during the storage, a partial hydrolysis of the active compound (adipoyl dichloride) can occur in the break-open phials containing organic phase. This hydrolysis is

accompanied by the release of hydrogen chloride, a gas which can cause an excess of pressure in the break-open phials containing organic phase.

In order to prevent any risk of eye projection during the opening of the break-open phials, it is imperative to wear safety goggles. We also advise you to keep this protection throughout the manipulation.

3.5 Collective protection

Since the reagents and the reaction products do not contain inflammable or toxic volatile solvents, no collective protection device is required during storage (safety cabinet, etc.) or use (hood, etc.).

3.6 Waste

The Nylon 6-6 synthesis kit does not contain any compounds which present a particular danger to the environment. However, for a responsible approach to environmental protection the products must be considered as laboratory chemical waste and treated as such. The quantities of reagents used in the kit have been reduced to minimize the volume of waste and the costs linked to its destruction.

4 Operating procedure

4.1 Principle of the synthesis

The nylon 6-6 preparation kit in break-open phial allows carrying the synthesis by interfacial polycondensation of a polyamide (in this case nylon 6-6) by making diamine (1,6-diaminohexane) in aqueous solution and acid dichloride (adipoyl dichloride) in organic solution molecules react between themselves.

The principle of the synthesis is based on the fact that the reaction takes place at the interface between the aqueous and organic phases that are immiscible. The reagents can only interact at the interface and react among themselves according to the following reaction (*Diagram 1*).


Diagram 1

The polycondensation reaction continues until the formation of a polyamide macromolecule, its formula is described below (*Diagram 2*).


Diagram 2

Diagram 2 shows an n degree nylon 6-6 macromolecule.

From a practical point of view the nylon is produced by catching the polyamide film formed in the interface and pulling it upwards while winding the film around a glass rod. The film is formed instantaneously as the nylon thread is wound and as long as there are reagent molecules in solution in the two phases.

4.2 Preparation of the experiment

At each work station assemble the recommended accessories (*paragraph 2*) and an break-open phial of each reagent.

A video is available on our site www.jeulin.fr.

Picture 2

4.2.1 Preparation of the reaction medium

Directions: to empty the bulbs, it is advisable maintaining the bulb as possible to the horizontal and practicing small pressure on the bulb.

Open the vial containing the aqueous phase carefully breaking the top (Picture 3). Pour the contents into a beaker of 50 or 100 mL. Repeat with the bulb containing the organic phase. Slowly pour the contents along the wall of the beaker (picture 4) while avoiding the eddies. The organic phase then floats on the aqueous phase.

Picture 3


Picture 4


4.2.2 Extraction of the nylon thread

Unfold the paper clip and form a hook at the end. Gently dip it in the reaction medium and catch the nylon film formed at the interface. **The quality of the final result depends on performing this step carefully.**

Then draw the film upwards over about ten centimeters. A nylon thread is formed (*picture 5*). Keep the thread in this position and bring the glass rod to mid-height of the thread. Then wind the five centimeters of the thread above the rod using the hook. Finally, wind the nylon thread by turning the glass rod on itself (*picture 6*).

Picture 5

Picture 6


4.2.3 Cleaning and drying of nylon thread

Once enough nylon has been extracted, cut the thread and dip the end of the glass rod into the second beaker previously filled with detergent solution (water with dishwashing liquid).

Shake the rod in the solution for a few minutes (*picture 7*) then dip it into the third beaker containing water. Place the piece of nylon on a sheet of absorbent paper and dry it by pressing on top with a spatula (*picture 8*).

Picture 7

Picture 8

5 After-Sales Service

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries