Mesures

Chronomètre

Measurement

Chronometer

Réf :
351 058

Français – p 1

English – p 5

Version : 7110

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2 Connections

The counter can be disconnected by exterior probes (mechanical contact or optoelectrical signal).

It always returns to open or close a contact between the **yellow** and **black** socket. This contact can be a physical mechanical switch, fork or optical barrier with an open collector output (JEULIN Timer ref. 353 026), or every optoelectrical probe with an open collector output.

For probe **1**, the sockets (yellow, black, red) situated on the left of the unit must be used, and for probe **2**, the sockets (yellow, black, red) situated on the right of the unit must be used.

For each optical probe the timer can supply a voltage of around 3.5 V with a current of 150 mA.

3 Working modes

There are 4 working modes:

3.1 Manual

This mode enables timing by manual release to be carried out at the start and end of the sequence (START / STOP button). It is possible to coagulate the display during the sequence (INTERMEDIARY TIME button), the timing continues and if this same button is pressed again the display is re-established with the time passes.

This mode is chosen by default when turned on.

The RESET button sets the display to zero when the timer is stopped.

3.2 2 Probes

This mode enables timing with an exterior release by 2 probes. Input 1 starts the timing and input 2 stops it. The RESET button stops the timer if it was started and sets it to zero.

Role of CLOSING / OPENING inverters:

Case of timer: In the OPENING mode, the timer is released when there is an opaque object placed in front of the optical beam, in the CLOSING mode, the timer is released when an opaque object is removed from the optical beam.

Case of **contact**: In the OPENING mode, the timer is released when a contact is opened between the yellow and black socket, in the CLOSING mode, the timer is released when a contact is closed between the yellow and black socket. Note that for a simple contact generating rebounds, the effect on the timing can appear to be the same.

3.3 1 Probe

This mode enables timing to be carried out with an exterior release by 1 probe adjusted to input 1. The first action on the probe releases the timer and the second action stops it.

The RESET button and CLOSING / OPENING inverter are the same functions, as when there are 2 PROBES.

The INTERMEDIARY TIME button has the effect of 'hiding' any movement in front of the probe when the timer is released (the timing continues). When the button is released, the next passing of an object in front of the probe stops the timing. Example of use of this mode: measurement of the pendulum period totalling several swings.

Note that due to the 'anti-rebound' detection, 2 successive releases must be separated by at least 150 milliseconds.

3.4 Timing 1

This mode enables the timing of an event from an optical fork or a contact to take place (on input 1).

Role of CLOSING / OPENING inverters:

Case of **timer**: In the OPENING mode, the incrementation of the timer is done so when there is an opaque object placed in front of the optical beam, in the CLOSING mode, the incrementation of the timer is done so when an opaque object is removed from the optical beam.

Case of **contact**: In the OPENING mode, the incrementation of the timer is done so when a contact is opened between the yellow and black socket, in the CLOSING mode, the incrementation of the timer is done so when a contact is closed between the yellow and black socket.

4 Maximum totalisation

In all modes, when the contents of the counter go above 20 000 points, 19 999s on the 20 s range or 199 99 s on the 200 s range, the display returns to zero and the 2 points flash. Press the RESET button to re-start the timer.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux