Méthode EXAO

ExAO method

Réf :
453 110

Français – p 1

English – p 4

Version : 7110

Module Angle SA
Angle Unit with analog output

1. Description

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
0 825 563 563*

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1. Description

The analog output angle unit, when attached to the electronic cone pulley (Part no.: 453109, not supplied), enables the computerization of angle measurements.

- 1 – 230 V / 12 V mains adaptor socket
- 2 – Push-button to select the rated capacities $-45^{\circ}/+45^{\circ}$ and $-180^{\circ}/+180^{\circ}$ with monitoring LED
- 3 – 0 calibration push-button for the 2 rated capacities
- 4 – Rated capacity calibration push-button $-45^{\circ}/+45^{\circ}$
- 5 – 6-pin DIN plug for connecting the pulley
- 6 – Analog output $-5\text{ V} / +5\text{ V}$
- 7 – Tilttable unbreakable support

2. Technical characteristics

Rated capacities: Selected by push-button (2):
 $-45^{\circ}/+45^{\circ}$, accuracy $\pm 2\%$
 linearity error: $\pm 1\%$
 $-180^{\circ}/+180^{\circ}$, accuracy $\pm 1\%$
 Resolution 0.75°

Calibration By push-button:
 $-45^{\circ}/+45^{\circ}$, at two points (3) and (4)
 $-180^{\circ}/+180^{\circ}$, at one point (3)

Dimensions (in mm) L. 100 x b. 100 x h. 40

3. Connection

The angle unit is powered by a 230 V /12 V mains adaptor (Part no. 281 512, not supplied). It is energised as soon as the adaptor is connected (1).

It processes signals coming from the cone pulley (Part no.: 453109, not supplied). The connection is made with a 6-pin DIN plug (5).

This device is equipped with an analog output (6) which helps the acquisition and processing of data from all interfaces and software compatible with this type of signals.

The characteristics of this analog output are:

- 5 V / + 5 V gives $-45^\circ / + 45^\circ$ or $-180^\circ / + 180^\circ$ according to the rated capacity selected,

Connection: safety banana sockets with 4 mm diameter.

4. Calibration

4.1 Heating time of the device

After powering on, the LED of the rated capacity – $45^\circ / +45^\circ$ comes on.

One must wait for about 1 min (the LED of the rated capacity flashes) during which the analog output is blocked at 0 V.

Once this time has elapsed, it is possible to make measurements. However, the device needs a heating time of about $\frac{1}{2}$ hour to obtain a stable zero value. If measurements are made during this period, it is recommended that the zero be calibrated before every acquisition.

This information is also valid when the rated capacity changes from 180° to 45° .

To change the rated capacity, press the corresponding push-button (2).

For the $-180^\circ / +180^\circ$ rated capacity there is no heating time. The measurements may be started directly.

4.2 Equipment calibration

Zero adjustment (whatever the rated capacity) :

Attach a pendulum on the pulley and wait till there is no movement. Press the zero adjustment push-button (3).

Angle adjustment (rated capacity – $45^\circ / + 45^\circ$):

Move the pendulum by $+ 45^\circ$ from its normal position (the screen printing on the pulley indicates the value of the angles).

Press the 45° adjustment push-button (4).

The device is calibrated and measurements may now be taken.

4.3 Computerized calibration

The “Science Studio” collection software (Généris 5 and 5 plus) help customize the analog output of the device in order to directly carry out the acquisitions of angles and not of voltage. Refer to the instruction manual of the software package (Acquisition chapter).

5. Experimentation

The angle unit, when attached to the electronic cone pulley (Part no.: 453 109) and to the simple and compound pendulum kit (Part no.: 222 022) helps perform many experiments on mechanical oscillations:

- Study of the natural period of a simple pendulum
- Study of the damping of a compound pendulum
- Verification of the isochronisms law for small oscillations of a simple pendulum
- Concept of amplitude, period, pseudo period, aperiodic and critical rates.

For further details regarding these uses, refer to the instructions of the simple and compound pendulum kit.

6. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux