Ex.A.O.

Tooxy[®] autonome

Data logging

Standalone Tooxy[®] console

Réf :
480 003

Français – p 1

English – p 6

Version : 5110

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0 825 563 563*

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1. Description

The standalone Tooxy® is a data logging interface capable of performing two measurements simultaneously making use of sensors with display.

It is designed to operate as a measurement device (read value directly on the sensor display) or as a data acquisition interface (connection to the computer through a USB cable).

It has rechargeable batteries so that it can be used without being connected to the mains power supply.

2. Use

2.1 Sensors

It is possible to connect the whole JEULIN range of sensors with display reference 482 XXX.

The sensor fits into one of the two ports. You will hear a "click" when the sensor is correctly connected.

2.2 Using the interface in standalone mode

Insert your chosen sensor in one of the ports.

A "click" informs you when the sensor is correctly inserted.

The measured value is displayed on the sensor.

2.3 Using the interface in connected to computer mode

Plug the USB cable into the interface.

Plug the USB cable into the computer.

The computer detects that new hardware is present. Wait until you see a message informing you that it is installed.

There is no need to reboot the computer.

2.4 Recharging

The standalone Tooxy® interface must be powered only using the power supply provided by Jeulin. (ref. 281 213).

Plug the transformer connector into the console.


Plug the power supply cable into the transformer.

Insert the plug into the mains power outlet.

The beginning of the charge is signaled by a light indicator flashing orange

The charge is finished when the indicating light is continuously orange.

Recommendations

You can use the standalone Tooxy® interface while the battery is being charged.

It is preferable to follow some recommendations to extend the life of your battery:

- Respect charge cycles:
 - o Recharge the battery only when it is almost completely discharged
 - o Do not disconnect and then reconnect the power supply several times during a charging cycle
- Recharge the interface before a long term storage period
- Recharge the interface at least once every 3 or 4 months, discharging it well before recharging it
- Store the device preferably in a location at a temperature between 10°C and 25°C

Information provided by the LED

Mains power supply	On/Off button	Battery state	Battery charge	LED state
Not connected	On	< Min charge	No	Red LED flashes for 30 seconds and then goes off
Not connected	On	> Min charge	No	Green
Connected	Off	< Min charge	Yes	Flashing orange
Connected	Off	> Min charge	Maintenance	Orange
Connected	On	< Min charge	Yes	Flashing orange

Complete recharge time : 3 hours

Endurance (full charge) : 7 hours

The endurance may vary depending on the consumption of the sensors used.

3. Technical characteristics

- Connection to the computer : USB
- 2 sensor inputs
- Max sampling rate : 10 kHz
- Resolution : 12 bits
- DC power supply 12 V, 1.25 A 

4. Precautions

The device must always be positioned facing the user. Do not block access to the device power supply connection.

The protection provided by the device may be compromised if it is used in a manner not specified in said manual.

Always disconnect the interface before performing any cleaning operation.

5. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux